

DO NOT MARK ON THIS QUESTION BOOKLET



Chemistry Lab Experiment Booklet



STOUT
UNIVERSITY OF WISCONSIN
WISCONSIN'S POLYTECHNIC UNIVERSITY

April 23rd, 2022

State Tournament

PART – 1 of 2

Instructions: This exam consists of a set of experiments based on oxidation-reduction reactions and aqueous solutions worth 100 points and a multiple choice exam worth 100 points (**200 points total**).

Students may NOT write on the Experiment or Exam Booklet. Students may **only mark on their answer sheets.**

Be strategic, figure out the way that you and your team can bank as many points as possible in the time given. Place the answers to the lab experiments and multiple-choice exam on the provided answer sheets. Answers not placed on the answer sheet will not be scored.

Ties will be broken by first the quality, accuracy, and completeness of the **experimental data and results**, followed by (if necessary) **selected multiple choice problems.**

You will need to perform the experiment and have your station completely cleaned up **and reset to their original positions** or a penalty of up to 10 % will be assessed against your score.


DO NOT MARK ON THIS QUESTION BOOKLET

Potentially Useful Information

$$1 \text{ mL} = 1 \text{ cm}^3$$

The density of water at room temperature, $\rho_{\text{H}_2\text{O}} = 1.00 \text{ g/mL}$

Table 1: Electrochemical Series		
Equilibrium (Oxidants \leftrightarrow Reductants)	E° (volts)	
Lithium: $\text{Li}^+(\text{aq}) + \text{e}^- \leftrightarrow \text{Li}(\text{s})$	-3.03	Metal Oxidizing Activity Increasing 
Potassium: $\text{K}^+(\text{aq}) + \text{e}^- \leftrightarrow \text{K}(\text{s})$	-2.92	
Calcium: $\text{Ca}^{2+}(\text{aq}) + 2\text{e}^- \leftrightarrow \text{Ca}(\text{s})$	-2.87	
Sodium: $\text{Na}^+(\text{aq}) + \text{e}^- \leftrightarrow \text{Na}(\text{s})$	-2.71	
Magnesium: $\text{Mg}^{2+}(\text{aq}) + 2\text{e}^- \leftrightarrow \text{Mg}(\text{s})$	-2.37	
Aluminum: $\text{Al}^{3+}(\text{aq}) + 3\text{e}^- \leftrightarrow \text{Al}(\text{s})$	-1.66	
Zinc: $\text{Zn}^{2+}(\text{aq}) + 2\text{e}^- \leftrightarrow \text{Zn}(\text{s})$	-0.76	
Iron: $\text{Fe}^{2+}(\text{aq}) + 2\text{e}^- \leftrightarrow \text{Fe}(\text{s})$	-0.44	
Lead: $\text{Pb}^{2+}(\text{aq}) + 2\text{e}^- \leftrightarrow \text{Pb}(\text{s})$	-0.13	
Hydrogen: $2\text{H}^+(\text{aq}) + 2\text{e}^- \leftrightarrow \text{H}_2(\text{g})$	0.00	
Copper: $\text{Cu}^{2+}(\text{aq}) + 2\text{e}^- \leftrightarrow \text{Cu}(\text{s})$	+0.34	
Silver: $\text{Ag}^+(\text{aq}) + \text{e}^- \leftrightarrow \text{Ag}(\text{s})$	+0.80	
Gold: $\text{Au}^{3+}(\text{aq}) + 3\text{e}^- \leftrightarrow \text{Au}(\text{s})$	+1.50	



DO NOT MARK ON THIS QUESTION BOOKLET

Potentially Useful Information

Solubility Table Common Ionic Compounds

	Group 1				Group 2			Transition Metals					
	NH ₄ ⁺	Li ⁺	Na ⁺	K ⁺	Mg ²⁺	Ca ²⁺	Ba ²⁺	Al ³⁺	Fe ³⁺	Cu ²⁺	Ag ⁺	Zn ²⁺	Pb ²⁺
F ⁻	sol	sol	sol	sol	insol	insol	sl sol	sol	sl sol	sol	sol	sol	insol
Cl ⁻	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol	insol	sol	sol
Br ⁻	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol	insol	sol	sol
I ⁻	sol	sol	sol	sol	sol	sol	sol	sol			insol	sol	insol
OH ⁻	sol	sol	sol	sol	insol	sl sol	sol	insol	insol	insol		insol	insol
S ²⁻	sol	sol	sol	sol		sol			insol	insol	insol	insol	insol
SO ₄ ²⁻	sol	sol	sol	sol	sol	sl sol	insol	sol	sol	sol	sl sol	sol	insol
CO ₃ ²⁻	sol	sol	sol	sol	insol	insol	insol				insol	insol	insol
NO ₃ ⁻	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol
PO ₄ ³⁻	sol	sol	sol	sol	insol	insol	insol	insol	insol	insol	insol	insol	insol
CrO ₄ ²⁻	sol	sol	sol	sol	sol	sol	insol		insol	insol	insol	insol	insol
CH ₃ CO ₂ ⁻	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol	sol


sol — soluble

sl sol — slightly soluble

insol — insoluble

(blank) — compound does not exist

FLINN
SCIENTIFIC, INC.
"Your Safer Source for Science Supplies"
© 2008 Flinn Scientific, Inc. All Rights Reserved.
AP6901



$$T(^{\circ}C) = [T(^{\circ}F) - 32^{\circ}F] \times \frac{5(^{\circ}C)}{9(^{\circ}F)}$$

$$T(^{\circ}F) = \left[T(^{\circ}C) \times \frac{9(^{\circ}F)}{5(^{\circ}C)} \right] + 32^{\circ}F$$

$$T(K) = [T(^{\circ}C) + 273.15(^{\circ}C)] \times \frac{1(K)}{1(^{\circ}C)}$$

$$T(^{\circ}C) = [T(K) - 273.15(K)] \times \frac{1(^{\circ}C)}{1(K)}$$

DO NOT MARK ON THIS QUESTION BOOKLET

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18																																																												
Hydrogen 1 H 1.008	Lithium 3 Li 6.94	Sodium 11 Na 22.990	Potassium 19 K 39.098	Rubidium 37 Rb 85.468	Cesium 55 Cs 132.91	Francium 87 Fr [223.02]	Beryllium 4 Be 9.0122	Magnesium 12 Mg 24.305	Calcium 20 Ca 40.078(4)	Strontium 38 Sr 87.62	Barium 56 Ba 137.33	Radium 88 Ra [226.03]	Titanium 22 Ti 47.887	Zirconium 40 Zr 91.224(2)	Hafnium 72 Hf 178.49(2)	Rutherfordium 104 Rf [261.10]	Scandium 21 Sc 44.956	Yttrium 39 Y 88.906	Lutetium 71 Lu 174.97	Lawrencium 103 Lr [260.10]	Vanadium 23 V 50.942	Niobium 41 Nb 92.906(2)	Tantalum 73 Ta 180.95	Dubnium 105 Db [261.10]	Chromium 24 Cr 51.996	Molybdenum 42 Mo 95.96(2)	Tungsten 74 W 183.84	Seaborgium 106 Sg [263.10]	Manganese 25 Mn 54.938	Technetium 43 Tc [98.906]	Rhenium 75 Re 186.21	Bohrium 107 Bh [264.10]	Cobalt 27 Co 58.933	Rhodium 45 Rh 101.07(2)	Osmium 76 Os 190.23(2)	Hassium 108 Hs [277.10]	Nickel 28 Ni 58.693	Palladium 46 Pd 106.42	Platinum 78 Pt 195.08	Darmstadtium 110 Ds [281.10]	Boron 5 B 10.81	Aluminum 13 Al 26.982	Gallium 31 Ga 69.723	Indium 49 In 114.82	Thallium 81 Tl 204.38	Ununtrium 113 Uut [284.18]	Carbon 6 C 12.011	Silicon 14 Si 28.085	Germanium 32 Ge 72.63	Tin 50 Sn 118.71	Lead 82 Pb 207.2	Flerovium 114 Fl [289.19]	Nitrogen 7 N 14.007	Phosphorus 15 P 30.974	Arsenic 33 As 74.922	Antimony 51 Sb 121.76	Bismuth 83 Bi 208.98	Ununpentium 115 Uup [288.10]	Oxygen 8 O 15.999	Sulfur 16 S 32.06	Selenium 34 Se 78.96(3)	Tellurium 52 Te 127.60(3)	Polonium 84 Po [209]	Livermorium 116 Lv [293]	Fluorine 9 F 18.998	Chlorine 17 Cl 35.45	Bromine 35 Br 79.904	Iodine 53 I 126.90	Astatine 85 At [210]	Ununseptium 117 Uus [294]	Helium 2 He 4.0026	Neon 10 Ne 20.180	Argon 18 Ar 39.948	Krypton 36 Kr 83.798(2)	Xenon 54 Xe 131.29	Radon 86 Rn [222.02]	Ununoctium 118 Uuo [294]


Key:
Element Name
Atomic number
Symbol
Atomic weight (mean relative mass)

Lanthanum 57 La 138.91	Cerium 58 Ce 140.12	Praseodymium 59 Pr 140.91	Neodymium 60 Nd 144.24	Promethium 61 Pm [144.91]	Samarium 62 Sm 150.36(2)	Europium 63 Eu 151.96	Gadolinium 64 Gd 157.25(3)	Terbium 65 Tb 158.93	Dysprosium 66 Dy 162.50	Hoiumium 67 Ho 164.93	Erbiumium 68 Er 167.26	Ytterbium 70 Yb 173.05	Actinium 89 Ac [227.03]	Thorium 90 Th 232.04	Protactinium 91 Pa 231.04	Uranium 92 U 238.03	Np 93 Np [237.05]	Plutonium 94 Pu [244.06]	Americium 95 Am [243.06]	Berkelium 97 Bk [247.07]	Californium 98 Cf [251.08]	Einsteinium 99 Es [252.08]	Fermiumium 100 Fm [257.10]	Mendelevium 101 Md [258.10]	Nobelium 102 No [259.10]	Ununseptium 117 Uus [294]	Ununseptium 117 Uus [294]	Ununoctium 118 Uuo [294]
--	-------------------------------------	---	--	---	--	---------------------------------------	--	--------------------------------------	---	---------------------------------------	--	--	---	--------------------------------------	---	-------------------------------------	-----------------------------------	--	--	--	--	--	--	---	--	---	---	--

*lanthanoids

**actinoids

Formula	Name
CATION: Positive Ion	
NH ₄ ⁺	ammonium ion
ANIONS: Negative Ions	
Based on a Group 4A element	
CN ⁻	cyanide ion
CH ₃ CO ₂ ⁻	acetate ion
CO ₃ ²⁻	carbonate ion
HCO ₃ ⁻	hydrogen carbonate ion (or bicarbonate ion)
Based on a Group 5A element	
NO ₂ ⁻	nitrite ion
NO ₃ ⁻	nitrate ion
PO ₄ ³⁻	phosphate ion
HPO ₄ ²⁻	hydrogen phosphate ion
H ₂ PO ₄ ⁻	dihydrogen phosphate ion
Based on a Group 6A element	
OH ⁻	hydroxide ion
SO ₃ ²⁻	sulfite ion
SO ₄ ²⁻	sulfate ion
H ₂ SO ₄ ⁻	hydrogen sulfate ion (or bisulfate ion)
Based on a Group 7A element	
ClO ⁻	hypochlorite ion
ClO ₂ ⁻	chlorite ion
ClO ₃ ⁻	chlorate ion
ClO ₄ ⁻	perchlorate ion
Based on a transition metal	
CrO ₄ ²⁻	chromate ion
Cr ₂ O ₇ ²⁻	dichromate ion
MnO ₄ ⁻	permanganate ion



DO NOT MARK ON THIS QUESTION BOOKLET

TO REDOX OR NOT TO REDOX, THAT IS THE QUESTION

Introduction

The purpose of this experiment is to determine the identity of an electrode paired with copper in an electrochemical cell and to determine if four additional reactions are RedOx reactions. All RedOx reactions will require full analysis. Several questions about the solution chemistry involved follow each experiment.

Materials:

12-well plate
Voltmeter
Squares of sandpaper
Copper electrode
1.0 M Copper (II) sulfate electrolyte (shared)
Unknown metal electrode
1.0 M Unknown metal sulfate electrolyte (shared)
Saturated potassium chloride (shared)
Tissue paper (for salt bridge)
Aluminum foil
Saturated sodium carbonate solution (shared)
1.0 M Hydrochloric acid (shared)
Magnesium filings
Plastic pipettes
Wood splint
Plastic tweezers
Deionized water (at sink)
Analytical balance
Pyrex beakers


Take careful notice of the layout of your station.

Points will be deducted (up to 10%) from your score if you do not clean and replace all of the items back to their original positions!


Safety Precautions:

1.0 M hydrochloric acid is a strong acid and very damaging to the eyes and skin! Wear chemical splash goggles to perform this experiment. Immediately rinse off any acid that gets on your skin with running water.

OPTIONAL: Gloves and a chemical-resistant apron or lab coat may be worn for this experiment as well if desired.

DO NOT MARK ON THIS QUESTION BOOKLET**Procedure 1****ELECTROCHEMICAL DETERMINATION OF ELECTRODE IDENTITY**

1. Locate your 12-well plate
2. Place 40 drops of 1.0 M **unknown metal sulfate** solution in the top left well (A1)
3. Place 40 drops of 1.0 M **copper (II) sulfate** solution to the left in well A2
4. Place 20 drops of **saturated potassium chloride** solution into the top right well (A4)
5. Prepare a twisted tissue paper salt bridge by cutting a 1 to 2 cm wide strip of tissue paper strip from a paper towel or kimwipe, twisting into a filament




6. Soak the paper filament in saturated potassium chloride solution, in well A4.
7. Obtain a copper and an UNKNOWN metal electrode. Sand your electrodes using the provided squares of sandpaper. Work over a sheet of paper to prevent scarring the benchtop. ALSO, be sure use the labeled “copper” and “unknown” sanding paper to avoid cross contamination of the electrodes with metal particles.




DO NOT MARK ON THIS QUESTION BOOKLET

8. Remove the tissue paper salt bridge from the potassium chloride solution, and drape between the left hand and right hand reservoirs of your electrochemical cell. Be sure the paper contacts the solution on both sides.



9. Place a copper strip in the RED alligator clip of the voltmeter
10. Place the unknown metal strip in the BLACK alligator clip of the voltmeter.
11. Use your voltmeter (set to **DC** volts, \bar{V}) to measure the voltage across the two electrodes, (notice: when the electrodes are attached properly, the voltage will be **POSITIVE**. If your voltage is negative, reverse the locations of the voltmeter leads. **BE CAREFUL NOT TO SPILL SOLUTIONS ONTO THE VOLTMETER!**)
12. Lower the copper into the copper electrolyte and the Unknown into the unknown metal sulfate electrolyte.




13. After the voltage stabilizes (may take several seconds, but usually less than a minute), record the **Cell Voltage** on your **Experimental Data Sheet**
14. The data at the very beginning of the experiment is more accurate than after letting the electrochemical cell run for awhile. Over time, the salt bridge causes the solutions to mix and can decrease the voltage.
15. Use the electrochemical series reference table at the beginning of this booklet to determine the identity of the UNKNOWN metal paired with copper in your electrochemical cell.
16. Your experimental voltage may be a bit different than the calculated theoretical voltage, but should be close enough to determine the likely identity of the UNKNOWN metal electrode.
17. Write the **Identity of the UNKNOWN Metal** electrode on your **Experimental Data Sheet**

DO NOT MARK ON THIS QUESTION BOOKLET

ANSWER THE FOLLOWING QUESTIONS ON YOUR Experimental Data Sheet

18. Write the **half-reaction** in the direction that it is occurring at the **COPPER** electrode:
19. Write the **half-reaction** in the direction that it is occurring at the now determined **UNKNOWN** electrode:
20. Write the **BALANCED NET Chemical Equation** representing the reaction occurring the electrochemical cell:
21. Which metal electrode is the **ANODE**?
22. Which metal electrode is the **CATHODE**?
23. Which metal **ATOM** (or ion) is being **OXIDIZED**?
24. Which metal **ATOM** (or ion) is being **REDUCED**?
25. Which metal **ATOM** (or ion) is the **OXIDIZING AGENT**?
26. Which metal **ATOM** (or ion) is the **REDUCING AGENT**?
27. **CONCEPT QUESTION:** How many **grams** of **copper (II) sulfate pentahydrate** would be required to prepare 100.0 mL of a 1.00 M solution of copper sulfate electrolyte?
28. **CONCEPT QUESTION:** What is the **PERCENT BY MASS** of **Cu²⁺ ion** in the resulting 1.00 M copper (II) sulfate solution? (The density of 1.00 M copper sulfate solution is 1.19 g/mL at 20°C)

**Procedure 2**

29. Procedure 2 will use the UNKNOWN electrolyte solution from your electrochemical cell solution from Procedure 1, the unknown electrolyte is the **sulfate salt** matching the electrode identity.
30. Disassemble your electrochemical cell (but do not dump your electrolyte solutions)
31. Remove the paper salt bridge and discard
32. Remove the electrodes and wipe them clean with a paper towel
33. Fill your 250 mL glass beaker about $\frac{3}{4}$ full with deionized water

DO NOT MARK ON THIS QUESTION BOOKLET

34. Using your plastic pipette, add **20 drops** of **deionized water** to well **B1**, directly below the (previously) UNKNOWN metal sulfate solution.
35. Place **5 drops** of **saturated sodium carbonate solution** into well **B1** and stir with your pipette
36. Rinse your pipette by repeatedly aspirating in your beaker of deionized water
37. Place several drops of the UNKNOWN metal sulfate electrolyte solution into well **B1** and stir with your pipette. You can mix by slowly aspirating the mixture several times with your plastic pipette.
38. Rinse your pipette by repeatedly aspirating in your beaker of deionized water
39. Make a note of your observations on your *Experimental Data Sheet*

ANSWER THE FOLLOWING QUESTIONS ON YOUR *Experimental Data Sheet*

40. Based on the identity of the UNKNOWN you determined in Procedure 1, propose and write a **BALANCED NET Chemical Equation** representing the overall reaction occurring in **Procedure 2**:
41. Is this reaction a **REDOX** reaction?
42. **If so**, write the **OXIDATION** and **REDUCTION half-reactions**. (**If not**, write the **TYPE** of reaction that **IS** occurring in both blanks of your answer sheet, and write “NA” for the answer to the following 4 questions.)
43. Which **ATOM** (or ion) is being **OXIDIZED**?
44. Which **ATOM** (or ion) is being **REDUCED**?
45. Which **ATOM** (or ion) is the **OXIDIZING AGENT**?
46. Which **ATOM** (or ion) is the **REDUCING AGENT**?
47. **CONCEPT QUESTION**: How many grams of the *now determined* UNKNOWN sulfate monohydrate would be required to prepare 100.0 mL of a 1.00 M solution of UNKNOWN sulfate electrolyte?

Procedure 3 – CAUTION: This procedure may create some fumes, avoid breathing with your face too close to the well plate . . .

48. Obtain a strip of aluminum foil and cut/tear off an approximately 1 cm x 1 cm square of foil
49. Place the foil into well **B2**, directly below the copper sulfate solution.




50. Using your plastic pipette, add **10 drops** of copper sulfate solution from well A2 to the aluminum foil in well **B2**

DO NOT MARK ON THIS QUESTION BOOKLET

51. Rinse your pipette by repeatedly aspirating in your beaker of deionized water
52. This reaction may take some time to complete.
53. Make a note of your observations on your *Experimental Data Sheet*

ANSWER THE FOLLOWING QUESTIONS ON YOUR *Experimental Data Sheet*

54. Look carefully at the contents of well B2, propose and write a **BALANCED NET Chemical Equation** representing the overall reaction occurring in **Procedure 3**:
55. Is this reaction a **REDOX** reaction?
56. **If so**, write the **OXIDATION** and **REDUCTION half-reactions**. (**If not**, write the **TYPE** of reaction that **IS** occurring in both blanks of your answer sheet, and write “NA” for the answer to the following 4 questions.)
57. Which **ATOM** (or ion) is being **OXIDIZED**?
58. Which **ATOM** (or ion) is being **REDUCED**?
59. Which **ATOM** (or ion) is the **OXIDIZING AGENT**?
60. Which **ATOM** (or ion) is the **REDUCING AGENT**?



Procedure 4 – CAUTION: *This procedure may create some fumes, avoid breathing with your face too close to the well plate . . .*


61. Using your plastic pipette, add **10 drops** of **deionized water** to well **A3**
62. Place **10 drops** of **1.0 M hydrochloric acid** into well **A3**.
63. Add **1 magnesium shaving** to well **A3**
64. Make a note of your observations on your *Experimental Data Sheet*
65. This reaction can take some time to complete.

ANSWER THE FOLLOWING QUESTIONS ON YOUR *Experimental Data Sheet*

66. Look carefully at the surface of the magnesium shavings, propose and write a **BALANCED NET Chemical Equation** representing the overall reaction occurring in **Procedure 4**:

DO NOT MARK ON THIS QUESTION BOOKLET

67. Is this reaction a **REDOX** reaction?
68. **If so**, write the **OXIDATION** and **REDUCTION half-reactions**. (**If not**, write the **TYPE** of reaction that **IS** occurring in both blanks of your answer sheet, and write “NA” for the answer to the following 4 questions.)
69. Which **ATOM** (or ion) is being **OXIDIZED**?
70. Which **ATOM** (or ion) is being **REDUCED**?
71. Which **ATOM** (or ion) is the **OXIDIZING AGENT**?
72. Which **ATOM** (or ion) is the **REDUCING AGENT**?
73. **CONCEPT QUESTION:** How many **mL** of 12.0 M concentrated hydrochloric acid would required to prepare 100.0 mL of 1.0 M hydrochloric acid in a 100 mL volumetric flask?

**Procedure 5**

74. Add **40 drops** of saturated sodium carbonate solution to well **B3**.
75. Watch carefully while adding **several drops** of **1.0 M hydrochloric acid** to the sodium carbonate solution.
76. Make a note of your observations on your *Experimental Data Sheet*

ANSWER THE FOLLOWING QUESTIONS ON YOUR Experimental Data Sheet

77. Propose and write a **BALANCED NET Chemical Equation** representing the overall reaction occurring in **Procedure 5**:
78. Is this reaction a **REDOX** reaction?
79. **If so**, write the **OXIDATION** and **REDUCTION half-reactions**. (**If not**, write the **TYPE** of reaction that **IS** occurring in both blanks of your answer sheet, and write “NA” for the answer to the following 4 questions.)
80. Which **ATOM** (or ion) is being **OXIDIZED**?
81. Which **ATOM** (or ion) is being **REDUCED**?
82. Which **ATOM** (or ion) is the **OXIDIZING AGENT**?
83. Which **ATOM** (or ion) is the **REDUCING AGENT**?

DO NOT MARK ON THIS QUESTION BOOKLET

84. **CONCEPT QUESTIONS:** Assuming the temperature in the lab is 70 °F, what is the approximate concentration of the saturated sodium carbonate in g/100g water? (see plot in “reference data” section at the beginning of the lab booklet)
85. What is the ***molality*** of the saturated sodium carbonate at this temperature?

Disposal:

- Use a pair of tweezers to remove any remaining aluminum foil and dispose of in trash can
- Place the copper and unknown metal electrodes in the Mixed Metal Waste beaker
- Use a pair of tweezers to put any remaining magnesium shaving in the Mixed Metal Waste Beaker.
- Carefully empty your well plate into the sink. All electrolytes can be disposed of down the drain.
- Rinse the well plate, be carefull not to splash water back up at you. This can happen if the water stream directly hits a well. Angle the well plate away from you down in the sink while rinsing
- Tap the rinsed well plate up side down onto several paper towels to remove any clinging droplets.
- **Thoroughly rinse all remaining beakers and pipettes, shake dry, and return ALL ITEMS to their ORIGINAL PLACES at your station.**

